

# *iScan 2 MULTI*

*As versatile as you are*




- For all large animal veterinarians
- Rugged and very light-weight
- Quick cleaning and disinfection
- Interchangeable probes



As versatile as you are


**iScan 2 MULTI**



Large, bright and easy to read display



Greater comfort of working with OLED goggles



Outstanding imaging quality

**Types of probes working with iScan 2 MULTI:**

**ABDOMINAL:**

- Convex 3.5 MHz R60
- Convex 5.0 MHz R50
- Linear 10.0 MHz L40
- Microconvex 6.5 MHz R15
- Microconvex 6.5 MHz R11

**RECTAL:**

- Linear Rectal 7.0 MHz L60
- Convex Rectal 5.0 MHz R60
- Microconvex Rectal 5.0 MHz R20

**SPECIALIST:**

- OPU 6.5 MHz R11
- Linear Backfat 3.7 MHz L150
- Linear Backfat 3.7 MHz L180
- Microconvex Endocavitar 6.5 MHz R10
- Linear Endocavitar 7.5 MHz L40



**Technical data**

|                                    |   |
|------------------------------------|---|
| <b>Presets</b>                     | Reproductive Tract, Pregnancy, Fetal sex, Late pregnancy, Custom preset to save                       |
| <b>Used technologies</b>           | LuciD – image noise reduction and contrast enhancement  |
| <b>Measurement</b>                 | Distance, Area, Volume, Grid, Age tables (Cow CRL, Cow BPD, DSG Horse, DO Horse, CRL Sheep, Lama BPD) |
| <b>Image refreshing</b>            | Up to 28 frames / sec   |
| <b>Memory</b>                      | 200 Image, 200 Cine loops   |
| <b>Data export</b>                 | Yes (to an external USB drive)  |
| <b>File format</b>                 | PNG, AVI  |
| <b>Display</b>                     | 7" LCD LED  |
| <b>Dimensions</b>                  | 256 x 165 x 69 mm   |
| <b>Device weight</b>               | 1520 g body only, 2640 g body with battery and probe  |
| <b>Probes</b>                      | Interchangeable electronic broadband  |
| <b>Recognition of probes</b>       | Automatic   |
| <b>Working time on one battery</b> | Up to 7 hours (continuous operation on a full charge)   |
| <b>Control</b>                     | Splash-proof keyboard   |
| <b>Accessories</b>                 | Shade cover, Goggles, Stand-off, Rectal probe extension for cattle or sheep.                          |



Interchangeable probes and battery



Rugged, light-weight, comfortable

